THE EUROPEAN DIRECTORATE FOR THE **QUALITY OF MEDICINES** & HEALTHCARE (EDQM)







Presentation of the CombiStats Online Application



Webinar 3 July 2024



Table of Content

Introduction

How to create my user account

What is the Workspace

CombiStats demo

Documentation

How to order CombiStats

Q&A session



The CombiStats Software – Desktop

- Statistical analysis of data of biological dilution assays in accordance with Ph. Eur. Chapter 5.3
- First release in 1999 to OMCLs, in 2005 to non-OMCLs
- Used by 394 organisations (960 licences) worldwide in 2023
- Users' expectations (top 3)
 - Audit trail
 - Login/password (user identification)
 - More interactivity (e.g. information/error messages)



CombiStats Online Application

- Go live: July 2024
 - Meet user's expectations, offering a new & interactive interface
 - Build a sustainable tool using a modern programming language



Technical requirements

- Computer with access to the Internet, Chrome or Edge browsers
- Acrobat Reader for viewing PDF reports, MS Excel for viewing extra analysis results



CombiStats Online Application

Flexible licence management system


- EDQM Store: 1 licence = n users/seats
- Possibility to assign and revoke users

Dedicated workspace for your organisation

- Access by login/password
- Folders & files creation, management, sharing
- Import of CombiStats Desktop 5.0 to 7.0 files

Modern & interactive interface

- Better flow of analysis options
- Enhanced graphical views
- Information and error messages
- Template, file protection
- PDF report, audit trail



Statistical Design & Analysis

Desktop options (wizard & advanced)

Size of	tab	les
---------	-----	-----

Number of tables, info. Lines

Number of pre-dilutions

Number of doses, replicates

Orientation of tables

Doses vertical/horizontal

Doses notation

Explicit/Implicit

Design

Completely randomised

Randomised block

Latin square > show design

ANOVA

No ANOVA to Complete

Equivalence

Model

Parallel lines, Slope ratio

3-parameter exponential curves

4- & 5-parameter sigmoid curves

Single-dose assay

Spearman-Kaerber approach

Transformation

Dropdown list/user-defined

Variance

Dropdown list/user-defined

Fixed parameters

slope, addition, multiplication

Prediction

Effective doses, inverse prediction

Weight function

Table of blank results (+ stats)

Multiple-dose (standard only) model

Subset analysis (PLA)

Confidence level: 80 to 99%

Plot of linearised data/reg. lines

Plot of studentised residuals

PDF report (audit trail)



Confidence levels, correlation/determination coefficients


	Desktop	Online application
Results	(fixed)	(range)
Common slope or intercept	90%	95% (80%-99%)
Potency estimates	95%	95% (80%-99%)
Effective doses, predictions	95%	95% (80%-99%)
Equivalence - individual slopes or intercepts	95%	95% (80%-99%)
Equivalence - differences or ratios	90%	90% (fixed)
Combination of results	95%	95% (80%-99%)


Desktop version (90% conf. level)

Common slope(factor): b = -111.255 (-115.612 to -106.898)

Correlation | r |: 0.991424

Online application (95% conf. level)





95% confidence level



How to create your user account


♣ Sign in


Your e-mail address	
Password	
Continue	>

Forgot your password?



Create your account







Your account has been created. You will shortly receive an e-mail with a link to define your password in order to finalise your registration.

You already have an EDQM account. Your access to this application has been granted. Please use your EDQM credentials to login.

Validate your account

Click the link received by email

[EDQM-Account] Creating your account - Please define your password -- [Compte EDQM] Création de votre compte - Veuillez définir votre mot de passe 🗾



From noreply-qual@edqm.eu on 2024-06-28 10:12

Details 1 Headers Plain text

Version française plus bas

*** This is an automatically generated e-mail, please do not reply ***

Dear My first name My las,

Welcome to the EDQM environment.

In order to activate your account, you will need to set your password. To do so, please follow the link below:

http://sign-qual.edqm.eu/validation/MzExOA/6mb-fc6404052bcf121c90cc/?

lang=en&redirect=aHR0cHM6Ly9jb21iaXN0YXRzLXF1YWwuZWRxbS5ldS9hY2NvdW50cy9sb2dpbg

Do not forward this link to other users, as it contains a token that is only valid for your e-mail address for 1 day.

Please be informed that your account will be deactivated if you do not set up your password and log in within 5 days.

For your information, your login is your e-mail address.

Permissions will be granted to your account, when necessary, for applications you need to access to. In this case, you will receive an automatic e-mail from the so-called applications.

Enter a password

Enter your password

Please enter your password twice so we can verify you typed it in correctly.

New password*

- Your password must contain at least 10 characters.
- Your password must contain at least 1 uppercase character.
- Your password must contain at least 1 digit.
- Your password must contain at least 1 special character among !@#\$%^&*(){}
- Your new password can not be identical to any of the 10 previously entered.
- Your password can't be a commonly used password.
- · Your password can't be entirely numeric.
- Your password can't be too similar to your other personal information.

New password confirmation*

Log in to finalise the activation

One last step

Password set. In order to finalise the activation of your account, you now need to log in.

Sign In




The Workspace

Secure & collaborative working environment




Workspace activation



Collaborative environment

Workspace (EDQM server)



Outside the workspace, exports of files





Records are exported as XML files (.epmx, .epax, and .epcx)



Analysis, combination reports as PDF.



Extra analysis results as Excel tables.

Actions

Create

Search, rename

Copy, move to

Export (zip file) / import

Add to my Favourites

Share (read-only or read and write)

Trash: items remain 30 days in the trash after deletion. Can be restored or deleted permanently



Ordering: CSID code for n seats (users)

From: Sales@edqm.eu <Sales@edqm.eu>

Subject: CombiStats - Order 000000

Your CombiStats licence - Your order YYYYYY - Our order confirmation OOOOOO

DO NOT REPLY TO THIS MESSAGE

Dear Customer,

Thank you for your order. Please find your licence activation key below:

COMBISTATS: STATISTICAL ANALYSIS OF DATA: access for N user(s)

CSID-AAAAA-BBBBBBBBBBBB-CCCCC-DDDD

The procedure for activating the licence key and registering users can be found here.

If applicable, please forward this message to the person responsible for managing CombiStats licences within your organisation.


If you experience any problems, do not hesitate to contact us via the **EDOM HelpDesk**.

Kind regards,

The EDQM Sales Section



Activation



Welcome to CombiStats

Web-based application for the Statistical Analysis of Biological Dilution Assay Results or Potency Assay Results

A licence is required to view content in this site. Click **here** to activate your licence.

Licence Activation

If you have a licence key, activate it here to create your workspace.

Otherwise, contact the licence manager within your organisation.

Licence key (CSID)*:

Once activated, your licence key will be valid for a year.

Activate



Users' access management

Activate

LICENCE

Workspace ID: 5BNDLNKJ

Number of access: 4 assigned - 3 available

Expiration date: 01 July, 2025

NEW LICENCE ACTIVATION

CSID code*:

In case of renewal, your licence key will be valid for a year. In case of extension of the number of seats, the expiration date remains unchanged.

Licence key	Access count	Subscription	Activated on	Activated by
CSID-5D0QG-XV5YK16Q4C-TPB57-H8P0	7	until 2025-07-01	01/07/2024	bao lalala

ASSIGNEES

E-mail*

Assign single user

☐ Can manage

Assign single

Choose file*

Assign multiple users in a click



Assign multiple

You can download a template to fill in here

E-mail	Assignee	Can manage workspace	
baothanh.nguyenvan@edqm.eu	bao lalala		Revoke



Demos

- Workspace
- Import of files created in Desktop version
 - Overview of analysis options & results
 - Subset analysis, show design
- CombiStats online application
 - Creation of an assay
 - Template, protection
 - PDF report, revision & audit trail
 - Combination of assay results



Useful links

Institutional website

https://www.edqm.eu/en/combistats-tm

User guide and FAQs

https://combistats.edgm.eu/user-manuals/combistats_user_guide.pdf/

https://combistats.edqm.eu/help/

• Privacy, security notices

https://combistats.edgm.eu/help/



Validation

Ensure continuity

Online application results = Desktop version results

Validation based on:

- all examples published in the **European Pharmacopoeia**, **chapter 5.3** 'Statistical Analysis of Results of Biological Assays and Tests'
- almost all numerical examples given by D.J. Finney in his standard work 'Statistical Method in Biological Assay, 3rd Edition, Grin, London (1978)'.

Calculated results and published results are equal

In addition, beta-version testing by several OMCLs in March 2024



Toward a smooth transition


- Go-live in July 2024
- Licence renewal in Oct. 2024 (webshop)
- CombiStats Desktop will be available until March 2025





How to order CombiStats

 From 5 July 2024, it will be possible to order CombiStats licence directly from EDQM webshop: https://store.edqm.eu/index.html



- Price: 600 euros per licence
- 2 types of licence:
 - Normal / renewable
 - Extension



The floor is yours!



Thank you for your attention



Stay connected with the EDQM

EDQM Newsletter: https://go.edqm.eu/Newsletter

LinkedIn: https://www.linkedin.com/company/edqm/

X: @edqm_news

Facebook: @EDQMCouncilofEurope

